

INSTAFLEX™ CONCENTRATE DUPLICATING MATERIAL



An Economical, All-Purpose, Fully Reversible Duplicating Hydrocolloid Concentrate

User friendly with a low melting temperature, minimal hold time, and easy cleanup. This universal duplicating material is ideal for producing both refractory and dental stone models, as well as molds for the fabrication of partial and full dentures (utilizing the fluid resin pouring technique).

1 gallon of Instaflex concentrate duplicating material makes 3 gallons. Just add water.

Typical Material Properties*

Melt Temperature	Idle Temperature	Hold Time	Pour Temperature	Color
195 °F (91 °C)	131 °F (55 °C)	20 minutes	136 °F (58 °C)	Pink

*These results are based on the testing methods, frequency and procedures of Ransom & Randolph or its approved suppliers. The levels referenced herein are only for general guidance and do not constitute a firm specification.

Application Instructions

1. Measure the required amount of Instaflex Concentrate duplicating material and tap water. **It is recommended that the full content of the Instaflex concentrate pail be used at one time.**
The ratio of Instaflex Concentrate duplicating material to water is:
 - 1 Part Instaflex Concentrate duplicating material : 1.66 Parts Water (by weight)
 - 1 Part Instaflex Concentrate duplicating material : 2 Parts Water (by volume)
2. Using a knife or scoop, cut the Instaflex Concentrate duplicating material into small pieces (approximately 1" cubes) and place into clean duplicating unit.
3. Pour warm tap water into duplicating unit.
4. Start duplicating unit cycle with melt temperature of 195 °F (91 °C) to fully melt the Instaflex Concentrate duplicating material.
5. Hold for 20 minutes.
6. Reduce unit temperature to pouring temperature of 136 °F (58 °C).

Reuse

To reuse Instaflex Concentrate duplicating material, thoroughly rinse used material under tap water, cut into small pieces (approximately 1" cubes), place into plastic bucket, and seal. When reusing Instaflex duplicating material, there is no need to add water. The water present after rinsing typically is sufficient replacement of any water lost in process.

Storage

Keep container closed when not in use. Store in a cool dry area.

Safety

Before handling this material, review the corresponding Safety Data Sheet (SDS).

Availability

Package Size	Color	Item Number
1 GAL (3.8 L)	Pink	021985

To order, visit shop.ransomrandolph.com or contact an authorized distributor. For a list of authorized distributors near you, visit



RANSOM & RANDOLPH

3535 Briarfield Boulevard | Maumee, OH 43537 USA
800.800.7496 | 419.865.9497 | 419.865.9997 (FAX)
www.ransom-randolph.com



Made in the USA

Ransom & Randolph GmbH
Leipziger Straße 40 | 04571 Rötha Germany
+49 342 06373999

Investing with Innovation™

INSTAFLEX™ CONCENTRATE DUPLICATING MATERIAL

ransom-randolph.com/distributors or contact your region's R&R Customer Service Team.

As the conditions or methods of use are beyond our control, we do not assume any responsibility and expressly disclaim any liability for any use of this product. Information contained herein is believed to be true and accurate but all statements or suggestions are made without warranty, expressed or implied, regarding accuracy of the information, the hazards connected with the use of the material or the results to be obtained from the use thereof. Compliance with all applicable federal, state, and local regulations remains the responsibility of the user. All potential liability related to the sale and use of this product is limited to the cost of the particular goods sold in their respective transactions.



RANSOM & RANDOLPH

3535 Briarfield Boulevard | Maumee, OH 43537 USA
800.800.7496 | 419.865.9497 | 419.865.9997 (FAX)
www.ransom-randolph.com



Made in the USA

Issue Date: July 26, 2024

Ransom & Randolph GmbH
Leipziger Straße 40 | 04571 Rötha Germany
+49 342 06373999

Investing with Innovation™