



A high-expanding, extra-strength investment ideal for casting printed patterns or traditional wax patterns using the flash-fire method.

Galaxy II (G2) universal investment is a phosphate-bonded dental investment designed for casting and ceramic pressing. It ensures consistent expansion, smooth surface quality, and accurate fit, making it suitable for various dental alloys and pressable ceramic systems.

Features

Phosphate-bonded investment for dental casting and ceramic pressing.

- Suitable for use with alloys and pressable ceramics.
- Supports consistent expansion for accurate fit.
- Helps produce smooth surface quality.
- Designed for dependable performance in demanding dental lab applications.
- Ideal for crown, bridge, and pressable restoration workflows.

Typical Material Properties*

Liquid:Powder Ratio

23:100

Working Time

5–6 minutes

**These results are based on the testing methods, frequency and procedures of Ransom & Randolph or its approved suppliers. The levels referenced herein are only for general guidance and do not constitute a firm specification.*

Availability

USA & INTERNATIONAL

Galaxy II investment

100 grams (100 count) kit.....1400-2
100 grams (100 count) box.....1400-2NL
15 lb. (6.8 kg) pail.....264151

Galaxy II investment liquid

16 fl. oz. (0.48 L) bottle.....1702
32 fl. oz. (0.95 L) bottle.....1402-11035

Each kit includes investment powder 100 (100g) pouches and (4) 16 fl. oz investment liquid.

To order, visit Shop.Ransom-Randolph.com, contact your region's R&R Customer Service Team, or contact an authorized distributor. For a list of authorized distributors, visit Ransom-Randolph.com/Distributors.

Rev. 2026-08-17